

```

no_elements = 1;
width=18.5e-3;
height=13e-3;
kerf = 0;
no_sub_x = 1;
no_sub_y = 1;
focus = [0 0 60]/1000;

Th = xdc_linear_array(no_elements, width, height, kerf, no_sub_x, no_sub_y, focus);
figure;
show_xdc_geir(Th, 1);

fs = 100e6; %sampling freq (100Mhz)
f0 = 2.5e6; % transducer center freq (2.5Mhz)
t0 = 1/f0;
dt = 1/fs; %sampling period
excitation = sin(2*pi*f0*(0:dt:1.5*t0));

figure;
plot(0:dt:1.5*t0, excitation);
xlabel("time (s)");
title("excitation pulse");

xdc_excitation(Th, excitation);

t_ir = -2/f0:1/fs:2/f0;
Bw = 0.6;
impulse_response = gauspuls(t_ir, f0, Bw);
figure;
plot(t_ir, impulse_response);
xlabel("time (s)");
title("impulse response");

freqz(impulse_response,1,1024,fs);

%sub problem2
%define a measurement point
x0=30e-3;
y0=10e-3;
z0=50e-3;
measure_point=[x0 y0 z0];
[h_x0, t_start]=calc_h(Th, measure_point);
figure;
plot(t_start+(0:length(h_x0)-1)*dt, h_x0);
xlabel("time (s)");
stitle = sprintf("spatial impulse response at (%2.2d %2.2d %2.2d)", measure_point);
title(stitle);

%sub problem3
no_elements = 1;
width=18.5e-3;
height=13e-3;
kerf = 0;
no_sub_x = 30;
no_sub_y = 30;
focus = [0 0 60]/1000;

Th = xdc_linear_array(no_elements, width, height, kerf, no_sub_x, no_sub_y, focus);
figure;

```

```

show_xdc_geir(Th, 1);

fs = 100e6; %sampling freq (100Mhz)
set_sampling(fs);
f0 = 2.5e6; % transducer center freq (2.5Mhz)
t0 = 1/f0;
dt = 1/fs; %sampling period
excitation = sin(2*pi*f0*(0:dt:1.5*t0));

figure;
plot(0:dt:1.5*t0, excitation);
xlabel("time (s)");
title("excitation pulse");

xdc_excitation(Th, excitation);

t_ir = -2/f0:1/fs:2/f0;
Bw = 0.6;
impulse_response = gauspuls(t_ir, f0, Bw);
figure;
plot(t_ir, impulse_response);
xlabel("time (s)");
title("impulse response");

freqz(impulse_response,1,1024,fs);

%
%define a measurement point
x0=30e-3;
y0=10e-3;
z0=50e-3;
measure_point=[x0 y0 z0];
[h_x0, t_start]=calc_h(Th, measure_point);
figure;
plot(t_start+(0:length(h_x0)-1)*dt, h_x0);
xlabel("time (s)");
stitle = sprintf("spatial impulse response at (%2.2d %2.2d %2.2d)", measure_point);
title(stitle);

```

Read rectangular data for plotting....

Plots aperture with physical element number...

Read rectangular data for plotting....

Plots aperture with physical element number...

Warning: Remember to set all pulses in apertures for the new sampling frequency







